IDENTSCARE LTDA - Av. Edgar Nelson Meister, 474 - Distrito Industrial - 89219-501 - Joinville - SC - Brasil / Braz. CNPJ / Tax ID: 05.106.945/0001-06 **INDÚSTRIA BRASILEIRA** / **BRAZILIAN INDUSTRY** - Responsável Técnico / Technical Authority: Friedrich Georg Mittelstädt CRC: 13100147 - SC - ANVISA nº 80172310079 - Marca / Brand: FGM® / Fabricado no Brasil / Made in Brazil

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EN Instructions for use

Glass-ceramic composite for CAD/CAM

For professional use only
Read all the information in this manual carefully before using the product. Keep it for consultation until the consumption of the product and when there is no more interaction with the last patient.

Product description
Brava is a composite cured at ideal conditions in the form of blocks for milling in automated lathes. The composite presents particles ranging from 40nm to 5.0µm, total inorganic load content in weight of 72% to 82% (52% to 60% in volume) in a highly cured resin matrix. The composite block has a conversion

degree of approximately 86%, significantly higher than the majority of the conventional composites (55% to 60%). This feature, along with the adequate distribution of particle sizes, results in a composite of high mechanical properties, mainly its hardness and resistance to wear, which contributes to the application in anterior and posterior teeth. Brava Block is able to provide excellent finishing with smooth and highly shiny surfaces.

Basic composition
Active ingredients: metacrylic monomers, initiators, co-initiators, stabilizers and silane.
Inactive ingredients: glass-ceramic load particles, silica and pigments.

Presentation
Kit: 5 units, 14L.
Low Translucency: A1-LT, A2-LT, A3-LT, A3,5-LT, B1-LT, C2-LT, Bleach.
High Translucency: A1-HT, A2-HT, A3-HT, B1-HT.

Product indication
Brava Block is recommended as permanent indirect restoration in the cases of:
· Total unitary crown on tooth
· Implant-supported crowns
· Inlays and Onlays
· Veneers

Precautions and contraindications

- For Odontology use only.
- In case of allergic reactions to the product discontinue its use and send the patient to a doctor.
- Manufacturing steps of the CAD/CAM systems are described in the instructions for use and in the manuals of each equipment manufacturer. The recommendations of the manufacturers should be followed.

- The milling process should be carried out through the irrigation strategy. Dry milling is not viable for Brava Block.
- Leave an internal space for the cement.
- The piece should be sand blasted on its internal surface.
- The piece should not be laser-etched.
- Do not use hydrogen peroxide solution, EDTA, bicarbonates or astringent products on the surface of the piece.
- Do not use temporary cements that contain eugenol.
- The use of self-adhesive cements is not recommended for the cementation of restorations using Brava Block due to the low interaction of those cements with the internal surface of the piece.
- Brava Block is supplied to the user in highly cured state.
- This material should not be fired under any circumstances during finishing or build up/Add.
- If near-pulp or pulp exposure has occurred, use a hydroxide layer first, followed by an application liner/base to minimize sensitivity.
- The material is not recommended for patients with clinical symptoms of bruxism.

Side effects
The product contains methacrylate monomers

that can cause irritation or mild allergic reactions in patients sensitive to these substances.

Instructions for use

- Secure the block in the machine following the instructions of the equipment manufacturer. Note: use the geometry that is compatible to the system to be used for the fitting of the block.
- Follow the phases described in the instructions for use of the different CAD/CAM systems for the milling process. Based on the digitalization data, model the restoration in the software for the creation of the data for the milling.
- Use the Brava Block strategy defined in the library of the equipment. When there is no Brava strategy, use the strategies for the lithium disilicate nano-hybrid composite or ceramic-hybrid.

Mechanical preparation of the tooth:
For the piece to have an adequate mechanical performance, it is necessary to follow certain orientations during the dental preparation, according to the procedure to be done:

Inlay and onlay:
The internal axial walls of the preparation should present expulsive angles of approximately 5° to allow correct seating of the pieces. Acute corners should be rounded. The occlusal reduction (cusps) should be of 1.5mm and the buccolingual distance should be at least 1.5mm. (Table 1 e Table 2).


Tabla 1: Inlay




Tabla 2: Onlay

Veneer:

The depth of the preparation on the facial surface varies in the different thirds of the teeth: in the cervical area, it is suggested to be 0.4mm, in the medium third the suggestion is 0.6mm. The incisal reduction recommended is from 0.5mm to 1.5mm. When "enveloping" the incisal border, all angles should be rounded (Table 3).


Tabla 3: Veneer

Total unitary crown:

The axial walls of the preparation must have approximately 5° of inclination and the wear depth must be of 1.0mm in the cervical area and 1.5mm in the medium third. The incisal/occlusal region should be reduced in 1.5mm. All of the preparation should present round angles (Table 4).

Tabla 4: Unitary crown

Occlusal adjustments in the piece and polishing:

Once ready, the piece should be fitted on the tooth before cementation. If any occlusal adjustments are needed, it can be done with fine grit burs. After those adjustments, polishing can be done with rubber or abrasive disks of medium and fine grit, felt disks Diamond (FGM) and the polishing paste Diamond Excel (FGM).

Adhesive preparation of the piece:

- Sandblast the internal part of the piece using aluminum oxide of up to 50 micrometers, with pressure of around 3 bar until the surface is matted.
- Proceed to the washing of the pieces in an ultrasonic cleaner for 180 seconds.
- Apply 70% alcohol and dry the surface for 30 seconds.
- Apply silane (preferably containing MDP) and let it react for 60 seconds. Remove excess with air blasts.
- Apply adhesive containing MDP (Ambar Universal APS, Ambar APS or Ambar - FGM) for 30 seconds and then apply air blast for 10 seconds to volatize the solvent and level the thickness of the adhesive layer. Then, light cure the whole surface for 20 seconds.

Adhesive preparation of the substrate:

On the tooth:

- Use of two-step adhesive (Ambar or Ambar APS - FGM):**
- Remove any dirt from the tooth preparation. If necessary, use Robinson brush and prophylactic paste for that purpose, followed by thorough washing with air/water blast. Isolate the tooth, preferably with absolute isolation;
 - Apply the 37% phosphoric acid etchant Condac 37

- (FGM) on enamel and dentin for 15 seconds. Then, wash it gently leaving the dentin humid;
- Apply the adhesive with friction all over the surface of the tooth, for 10 seconds, in two layers. Then, apply air blast gently to spread the adhesive better and reduce the adhesive layer thickness. Light cure for 20seconds.

. Use of the one-step Ambar Universal APS (FGM) universal adhesive, compatible with the FGM line of resinous cements:

- This adhesive may be used associated or not to the acid etching of the enamel and dentin. If the professional chooses to etch, the protocol is identical to the one described above for Ambar or Ambar APS (FGM) adhesives. In case the professional chooses NOT to apply the acid etching, the protocol to be followed is here:
- Remove any dirt from the tooth preparation. If necessary, use Robinson brush and prophylactic paste followed by thorough washing with air/water blast. Isolate the tooth, preferably with absolute isolation;
 - Apply the adhesive with friction all over the surface of the tooth, for 10 seconds, in two layers. Then, apply air gently to spread the adhesive better and reduce the adhesive layer thickness. Light cure for 20s.

On metallic surfaces (prosthetic components over dental implants):

- Apply a "primer" for metal (preferably with MDP) on the prosthetic component already activated/installed on the dental implant.

Note: it is highly recommended that the excess cement is precisely emptied with an analog of the prosthetic component (extra-orally), immediately before the piece is placed in the mouth on the prosthetic component.

Cementing the piece to the tooth:

- With the prosthetic piece and the tooth prepared, in an isolated operatory field free from fluids, the piece is cemented. For that, manipulate the cement or use the dual-body syringe of Allcem or Allcem Core (FGM) according to the instructions for use of the chosen product;
- Apply the cement inside the prosthetic piece and place it on the tooth. The excess cement that may overflow must be removed. Light cure each surface of the tooth (labial, palatal/lingual and occlusal) for 20 seconds each, totalizing 60 seconds of light application. After that, wait 7 to 10 minutes for the cement to cure chemically.
- Check for proximal excess with dental floss and, if necessary, remove remaining resin/adhesive excess.

Repairing the piece after cementation:

- If, at any time, the professional needs to repair a prosthetic piece due to fractures or chipping, the following is suggested:
- Remove any dirt from the surface of the restoration area. If necessary, use Robinson brush and prophylactic paste. Isolate the tooth, preferably with total dam isolation.
 - Roughen the area to be repaired with diamond tipped drill bits of any granulation, wash and dry the surface;
 - Apply silane and adhesive with MDP on the piece, according to the product instructions;
 - Do the necessary repair using conventional composite resin (Llis, Opallis or Vittra APS) and light cure the increments according to the product instructions;
 - Proceed with the finishing and polishing according to

the steps of a conventional restoration.

Conservation and storage

Protect from direct sunlight. Store the product between the temperatures of 5°C and 35°C / 41°F and 95°F.

Warnings

- Make sure to use water during the milling process (wet operation). Dry operation is not suitable for Brava Block.
- The product is for professional use only.
- Do not use the product after its expiration date.
- Follow your country's legislation when disposing it.
- Keep out of the reach of children.
- Avoid the use in patients with a history of allergy to the composite resins containing methacrylates. In case of allergy, suspend the use and send the patient to a doctor.

This material was manufactured for dental use only and must be manipulated according to the instructions for use. The manufacturer is not responsible for damages caused for other use or incorrect manipulation. Additionally, the user is obligated to prove, before the use and under his/her responsibility, if this material is compatible with the desirable utilization, mainly in cases in which such utilization is not indicated in this manual. Data description do not constitute any type of warranty and, for this reason, are not binding.

Legenda de símbolos Meaning of symbols Leyenda de los símbolos	
	Data de fabricação Date of manufacture Fecha de fabricación
	Número do lote Batch Number Número del lote
	Prazo de validade Shelf life Fecha de validez
	Limites de temperatura Temperature limits Limites de temperatura
	Não utilizar se a embalagem estiver danificada Do not use if packaging is damaged No lo utilice si el empaque está dañado
	Conservar seco Keep it dry Consérvelo seco
	Mantter afastado da luz solar Keep protected from sunlight Manténgalo lejos de la luz solar
	Fabricante Manufacturer Fabricante
	Consulte as instruções de utilização Refer to instructions for use Consulte las instrucciones de utilización

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