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Characterization System Frantins

ACRYLIC RESIN

Indicated for the characterization of prosthesis according to the patient's gingiva staining.

Acrylic Resin Termopolimerizável, developed for works of characterization of gums through the method and Characterization Frantins. It is a work system because it contains a written method that allows for the repetition of the characterizations in a scientific way.

The method allows the customization of the prosthesis according to the patient's gingiva staining.

Type 1; class 1

Composition:

Powder: Polymethyl, benzoyl peroxide and biocompatible pigments.

Liquids: NO Crosslink: Methyl methacrylate and inhibitor.

with Crosslink :: Liquid Methyl methacrylate, EDMA (Crosslink) and inhibitor.

Precautions:

- * Read the INSTRUCTIONS on the label before handling .
- * For dental use only.
- * Flammable Liquid, non-toxic (Flash Point: 10°C / 50°).
- * Mix Air / powder (resin) becomes flammable.
- * It is recommended handling in a ventilated area wearing Nitrile gloves and protection glasses.
- * Pregnant women, Infants and sensitive people should NOT handle it.
- * Keep away from light and heat.
- * Do not allow polymerized material get in contact with the skin of the patient. The heat Polymerization CAN cause burns.

* **Colors / Powder: No veins: red rose, black rose, light pink, brown rose, purple rose, medium rose**
Veins Color: Red and Purple.

* **Packaging: full box- Powder: 06 Bottles of 15g,.**

Veins: 02 bottles of 01g

Individual kit box- Powder: 08 bottles of 15g with single color or 08 bottles with veins.

SHAKE THE BOTTLE BEFORE USING

*** Preparation:**

Before starting work, the parts to be used must be completely clean and free of traces of wax and fat, keep clean hands so there is no color change.

RECOMMENDATIONS: For every each 14grs. dose(19ml) of resin use 6,5 ml. of liquid. Tests must be done between 15°C and 30°C. Use only TERMODURAL liquid with crosslink. After the dosage of resin and liquid TERMODURAL and the mixture well homogeneous, the plastic fase begins(when the resin comes off the container wall)10 min.after the dosage 25°C mold the resin with the plaster already isolated. Press carefully until 1000kgf.(1ton.) for 1 minute. Remove the excess and press again until 1200kgf.(1.2ton.) and let it rest, pressed, for 20min before polymerizer.

Important; the ideal to increase translucency and porosity would be to keep for 2 hours before polymerizer.

***Conventional Polymerization:** On a pan with 3 liters of cold water, place the muffle, turn the stove on and let it reach the temperature of 70°C. Keep this temperature for 25 to 30 min. and after increase the temperature and let it reach 100°C (boiling point) and keep it for 1:30h. Let the muffle to cool in the water (room temperature) for 20 min. before start the deflasking.

Thermo-Pneumatic Polymerization: Place the muffle in the polymerizer, covering with water. Close and inject compressed air (60 pounds). Turn on the electric power until it reaches 120°C. Turn it off when the temperature is reached and let it cool down until 90°C. Once this temperature is reached, turn it on again until it reaches 120°C. Turn it off and let it cool down until 60°C. Make sure there is no pressure inside the polymerizer, open it and withdraw the muffle to begin the deflasking. Be careful, the muffle must be cold so it can not damage teeth and prosthesis.

Finishing and Polishing:

Finishing is started with bits and mills, passes to abrasive stones and it is concluded with fine texture sand papers.

Note: Thermodural acrylic liquid is presented with two options, normal and Crosslink to which allows great impermeability, preventing impregnations, unpleasant odors and favoring a resistant prosthesis, unalterable and without irritation in the oral environment.

1. Composition and information about ingredients.

Ingredients	CAS Number	Risk Symbols	Concentration
Polymethylmethacrylate	9011-14-7	N/A	N/A

2. Technical Responsible.

Regiane do Nascimento Rodrigues – CRQ 04163305
Anvisa Registration: Class: II 80984220002

3. Other informations.

N/A Not Applicable N/C Not Classified N/D Not Available

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