TERMODURAL

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RESINA ACRÍLICA

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Acrylic Resin thermo-polymerizable(powder/liquid): for dental prosthesis, applicable for Prosthesis Total Removable, Prosthesis Partial Removable and Prosthesis Total Fixed.

Type 1; Class 1

Composition:

Powder: Polymethacrylate, Benzoyl Peroxide, Biocompatible Pigments.

Liquid: No Crosslink: Methylmethacrylate

Liquid: With Crosslink: Methylmetacrylate, EDMA(crosslink) and inhibitor.

Precautions:

- Read the handling instructions described on the label before using it.
- For dental use only
- Flammable liquid, non-toxic material (Flash Point: 10°C / 50°C).
- Air/Powder mixture (Resin) can form flammable.
- Handling the product is recommended in ventilated areas, wearing nitrile gloves and protective goggles.
- Pregnant women, nursing mothers and people with known allergies should not handle this product.
- Keep protected from heat and light.
- Do not let a NON polymerized material get on the patient's skin. The heat from polymerization can cause burns.

PRESENTATION

PACKAGING:

POWDER:

Plastic bottles – 80g, 225g, 450g, 1000g

Paper barrels – 2.250g, 5.000g, 10.000g, 25.000g, 50.000g

LIQUID:

No Crosslink: glass bottles: 60ml, 120ml, 250ml, 500ml, 1000ml and 5 liters (gallon).

Crosslink: glass bottles: 60ml, 120ml, 250ml, 500ml, 1000ml and 5 liters (gallon).

COLORS/POWDER:

No veins: Colorless, Cristal and Palate

With veins: Light pink, Medium Pink, Medium-B, Dark Pink, Black and Special.

SHAKE WELL BEFORE USE.

•PREPARATION:

RECOMMENDATIONS: For every each 14grs. dose(19ml) of resin use 6,5 ml. of liquid. Tests must be done between 15°C and 30°C. Use only TERMODURAL liquid with crosslink.

After the dosage of resin and liquid TERMODURAL and the mixture well homogeneous, the plastic fase begins(when the resin comes off the container wall)10 min.after the dosage 25°C mold the resin with the plaster already isolated. Press carefully until 1000kgf.(1ton.) for 1 minute. Remove the excess and press again until 1200kgf.(1.2ton.) and let it rest, pressed, for 20min. before polymerizer.

•**Conventional Polymerization:** On a pan with 3 liters of cold water, place the muffle, turn the stove on and let it reach the temperature of 70°C. Keep this temperature for 25 to 30 min. and after increase the temperature and let it reach 100°C(boiling point)and keep it for 1:30h. Let the muffle to cool in the water(room temperature) for 20 min. before start deflasking.

Thermo-Pneumatic Polymerization: Place the muffle in the polymerizer, covering with water. Close and inject compressed air(60 pounds). Turn on the electric power until it reaches 120°C. Turn it off when the temperature is reached and let it cool down until 90°C. Once this temperature is reached, turn it on again until it reaches 120°C. Turn it off and let it coll down until 60°C. Make sure there is no pressure inside the polymerizer, open it and withdraw the muffle to begin the deflasking.

- Solubility in water, insoluble.

1.Composition and Informations about the raw materials.

Ingredients	CAS Number	Risk Symbols	Concentratrion
Polymethylmethacrylate	9011-14-7	N/A	N/A

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•**Anvisa Registration:** Class: II 80984220002

3.Other informations.

N/A Not Aplicable N/C Not Classified N/D Not Disponíble.

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