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Self-cured acrylic resin (powder / liquid) for repair, relining and additions. **Tipe 2; Class 1.**

Composition:

Powder: Polymethylmethacrylate, Benzoyl Peroxid, Biocompatible pigments

Liquid: Methylmethacrylate, DMT, Inhibitor

***Colors/ Powder:** colorless and pink.

Packaging

***Packaging Powder:** Polyethylene Bottles: 80g, 225g, 450g, 1,000g and 5,000g

Aluminized Refill Sachet: 80g, 225g, 450g, 1,000g and 5,000g

Liquid:

Glass: 60 ml, 120 ml, 250 ml, 500 ml, 1000 ml

Refill Bottles: 60 ml, 120 ml, 250 ml, 500 ml, 1000 ml and 5 liters (blue gallon).

Precautions:

*Read the instructions before handling.

*For dental use only.

* Fluid flammable, non-toxic material (Flash point: 10°C / 50°F)

*Mixture air/powder can become flammable

* Handling of the product is recommended in ventilated places, using nitrile gloves and protective goggles.

* Handling is not recommended for pregnant women, nursing mothers and sensitive people to this type of product.

*Keep away from sun light and heat.

*Do not let material NO POLYMERIZED get in contact with skin, the heat from polymerization can cause burns.

Instructions:

Shake the bottle (liquid) before using it.

*Preparation:

Before starting any work, the parts to be used should be completely clean and free from traces of wax, dust or fat, so that the chemical union between the resins may be perfect.

Recommendations: The ideal for each dose of resin (14g) uses 6,5 ml of liquid. Tests must be realized at temperatures between 15°C e 30°C. Use FAST (self) liquid.

Place the necessary amount of liquid (by the dripping method) in a jar of glass or silicon according with the amount of acrylic. Remove any excess. Dispose it and homogenize carefully, avoiding the incorporation of air bubbles. Observe so that the mixture is fluid, however consistent (do not touch on the prepared surface). Cover with polyethylene film (do not use Cellophane). Do not use the resin if stops sliding from the spatula (between 5 and 6 minutes after dosage at 25°C). Press with 1000kgf (01ton.) for 1 minute.

Remove the excess then the film, close with manual clamp.

Polymerization:

Take the resin to the polymerizer with enough water at room temperature. Apply pressure (20lbs.) of compressed air for 25 minutes. Wait for total polymerization to begin the polishment.

Finishing and Polishing

Finishing is started with bits and mills, passes to abrasive stones and it is concluded with fine texture sand papers.

1.Composition and Informations about the raw materials.

Ingredients	CAS Number	Risk Symbols	Concentration
Polymethylmethacrylate	9011-14-7	N/A	N/A

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