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FRANTINS INDUSTRY AND TRADE OF DENTAL PRODUCTS LTDA-  
EPP  
BRAZIL – SÃO PAULO – PIRASSUNUNGA – CURITIBA STREET, 4721  
– V. BELMIRO/ ZIP CODE 13633-515  
**PHONE: (19) 3097-1209**

## INDICATIONS:

**ACRYLIC RESIN THERMO-POLYMERIZABLE SPECIFICALLY FOR  
WORKS USING MICRO WAVES ENERGY.**

Type1; Class 1

Resin of high molecular weight, sustain high temperatures, for this technique **only MICRON liquid** must be used, where all the advantages are concentrated for conventional jobs.

## COMPOSITION:

**Powder:** Polymethyl metacrylate, Benzoyl Peroxide and Biocompatible Pigments.

**Liquid:** No Crosslink: Methylmetacrilate and Inhibitor.

**Liquid:** With Crosslink: Methylmetacrilate, EDMA(crosslink) and Inhibitor.

## Precautions:

- Read the handling instructions described on the label before using it.
- For dental use only
- Flammable liquid, non-toxic material (Flash Point: 10°C / 50°C).
- Air/Powder mixture (Resin) can form flammable.
- Handling the product is recommended in ventilated areas, wearing nitrile gloves and protective goggles.
- Pregnant women, nursing mothers and people with known allergies should not handle this product.
- Keep protected from heat and light.
- Do not let a NON polymerized material get on the patient's skin. The heat from polymerization can cause burns.

## Packaging:

Resin Micron plastic bottles: 225g and 450g

Liquid Micron with crosslink, glass bottles: 120 ml and 250 ml.

## SHAKE THE POWDER BEFORE USING IT.

### .Preparation:

Utilize ISOTINS insulator to guarantee prosthesis free of plaster incrustations, making deflasking and polymerization of the acrylic resin easier. Hands and containers must be clean to avoid color alteration. The model in wax is prepared and placed in the muffle normally. The wax is melted and the muffle cavity washed with hot water and detergent or according by the professional method. The cavity and the plaster base are insulated by an Alginate based insulating agent.

**Dosing and Plasticization:**

Place 6.5ml of the monomer in the mixing container, usually glass. Add 14g (19ml) of the resin in the desired color and with a plastic spatula mix the resin and the monomer until the homogenization is perfect. Wait for the plastic phase, which occurs when the material is set loose from the container's walls.

**NOTE:** Keep the muffle humid, especially if the work will be finished in the next day.

**Pressing:**

In the plastic phase, the resin must be placed in the muffle insulated in the preferred way. The resin is covered by a plastic film (high density polyethylene, do not use cellophane). Take to the muffle to press and perform pressing tests (slowly up to 500/1000kgf). After removing the film and the excess material, close the muffle and perform a definitive pressing (at the most 1000kgf). For prosthesis for at most 10mm or for prosthesis done with the use of silicon wall, wait for 15 minutes in the press or outside of it, before starting the polymerization in the microwave. For prosthesis above 10mm of thickness, the waiting time recommended is 2 hours. For prosthesis above 20 mm, 2 hours waiting is indispensable.

**Alternative:**

Instead of waiting for the plastic phase, some professionals pour a liquefied resin in the muffle, cover it with plastic film and wait for the plastic phase to come inside of the muffle before initiate the pressing process.

**Observation:**

Polymerization cycles by Frantins Methods of Characterization can also be used, according to the product manual instructions.

**Deflasking:**

Deflasking is done normally, make sure the muffle is cold and be careful not to damage the teeth or the prosthesis.

# MICRON

ACRYLIC RESIN

## Technical Datasheet English

### Polimerization:

Polymerization must be done extremely right, obeying the parameters by the table below:

| To 500w oven                 |     |                       |
|------------------------------|-----|-----------------------|
| INITIAL<br>STAGE             | 20  | With 20/30% of power  |
|                              | MIN | or Medium             |
| Low/defrost                  |     |                       |
| FINAL<br>STAGE               | +5  | With 80/100% of Power |
|                              | MIN | Or médium High -High  |
| To Owen from 800 w to 1100 w |     |                       |
| INITIAL<br>ISTAGE            | 20  | With 10/20% of power  |
|                              | MIN | Or low/defrost        |
| FINAL<br>STAGE               | +5  | With 50/60% of Power  |
|                              | MIN | or Medium             |
| To owen from 1200 w to 1400w |     |                       |
| INITIAL<br>STAGE             | 20  | 10% Power             |
|                              | MIN | or Low/Defrost        |
| FINAL<br>STAGE               | +5  | With 30/40% of Power  |
|                              | MIN | or Medium Low-Medium  |

### 1.Composition and Informations about the raw materials.

| Ingredients            | CAS Number | Risk Symbols | Concentration |
|------------------------|------------|--------------|---------------|
| Polymethylmethacrylate | 9011-14-7  | N/A          | N/A           |

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